

TechNiche LLC 2575 Pioneer Ave, Suite 101 Vista, California 92081 888.823.2665 techniche-intl.com

## Product #6865

TechNiche Phase Change Cooling H20 COOLING INSERTS

- → Activated by placing in the freezer
- → Use additional inserts for continuous relief
- → Water based Cooling Inserts maintain a temperature of 32°F (0°C) for up to one hour
- → Includes 4 Water Based Cooling Inserts per set → Weight: 6 lbs
- → Compatible with: 6826, 6626, 6626HV

H20 Cooling inerts uses a water based liquid that freezes to a temperature of  $32^{\circ}F$  (0°C)and remains at that temperature for up to 1 hour.

## HOW TO USE :

- 1. Simply place H20 inserts in freezer for 45 minutes
- 2. Insert into garment pockets
- 3. Garment is ready to wear



## FOLLOW US ON SOCIAL MEDIA



To Stay Up To Date On The Latest Deals & News!